



Reusable strap tag attachment

Please follow these instructions to ensure optimum performance from your ChronoTrack timing tag.



If you have difficulty following these instructions please contact a Timing Official for assistance

 Check that the number on the label is your race number. The tag is programmed with the number on the label.



Insert velcro strap into right hand slot on your tag



slide tag and thread strap through left hand slot



Slide tag along strap and pull strap firm ...



then firmly press velcro closed.



Wrap strap around ankle and **firmly press closed**.



The tag in place ready to record your times.



Please note **your ChronoTrack timing tag must be returned after the event.**



Failure to return your tag will result in you being charged for a replacement tag.



If you lose your tag during the event please advise an Event Official to enable a manual time to be recorded.